
# REFLOW MEDICAL

## WINGMAN CROSSING CATHETER Instructions for Use

Caution: This device is restricted to use by or on the order of a physician.

### Description

The Wingman Crossing Catheter consists of a support catheter, with an extendable tip, and activating handle. The crossing catheter is tracked over a guidewire to the occluded lesion. The extendable tip is in the retracted position during tracking over the guidewire. Upon reaching the occlusion, the tip is advanced forward, rotated, and retracted using the advancer on the handle. The advancement and retraction of the tip supports and allows the wire to penetrate and cross the lesion. This step can be repeated until the occlusion is crossed. When the tip is in the normal mode it is retracted just flush to the support catheter. When it is activated it is extended just distal to the catheter body. The crossing catheter will also allow for exchange of guidewires and other interventional devices, and provide a conduit for delivery of saline solutions or diagnostic contrast.

### Indications for Use

The Wingman Crossing Catheter is intended to be used in conjunction with steerable guidewires to access discrete regions of the peripheral vasculature. It may be used to facilitate placement and exchange of guidewires and other interventional devices, and provide a conduit for delivery of saline solutions or diagnostic contrast.

### Contraindications

The Wingman Crossing Catheter is contraindicated for use in the coronary and cerebral vasculature.

### Warnings

Hydrophilic wires prone to excessive swelling (e.g. ZipWire) should not be used with the Wingman Crossing Catheter.

Single Use only. Do not reuse/sterilize. Reusing the device could result in compromised device performance, cross-infection and other safety related hazards.

Do not use if device is open or packaging is damaged

Never advance, withdraw or rotate an intravascular device against resistance until the cause is determined by fluoroscopy.

### Directions for Use

Note: Follow instructions for use on all equipment to be used with the Wingman Crossing Catheter

- Preparation for Use:
  - Using sterile technique, remove the Crossing Catheter dispenser coil from its packaging and transfer it to the sterile field.
  - Remove the Crossing Catheter from the dispenser coil and inspect for any bends or kinks.
  - Attach a 10ml syringe filled with sterile, heparinized saline to the luer-lock guidewire entry port of the Crossing Catheter and thoroughly flush the catheter.
  - Prior to use, wet the distal 40cm of the Crossing Catheter with Heparinized saline solution to activate the hydrophilic coating.
- Insertion
  - Through a previously inserted, appropriately sized guiding catheter or introducer sheath, introduce the distal end of the Crossing Catheter over an appropriately sized guidewire (see specifications) using standard technique.
- Advancement
  - Use fluoroscopic guidance when advancing the Crossing Catheter to the desired location within the vasculature.
- Activation
  - Continue advancing the Crossing Catheter until lesion is reached. Using fluoroscopic guidance, confirm that Crossing Catheter and guidewire are at the lesion.
  - Once position is confirmed, slightly retract the guidewire within the CROSSING catheter and advance the extendable tip by depressing the handle down while turning - 90 degrees in a clockwise motion. This will extend the tip beyond the catheter, penetrating the lesion and providing support for the guidewire.  
**CAUTION: DO NOT turn the handle more than 90 degrees during advancement of the extendable tip. The tip will not advance further by turning more than 90 degrees; continuing to rotate the handle beyond 90 degrees may lead to device breakage.**  
Note: The handle should only be directed in a forward and backward motion. A side to side motion will not actuate the tip. Note: Excessive grip should not be placed on the outer shaft during advancement, as this will restrict the movement of the extendable tip.  
**CAUTION: Never advance, withdraw or rotate the CROSSING catheter against resistance until the cause is determined by fluoroscopy.**
  - Once penetrated, advance the guidewire through the lesion and gently retract the Crossing Catheter.
  - Repeat as needed until the Crossing Catheter and guidewire have passed the lesion.
- Infusion
  - To perform infusion, withdraw the guidewire and reference the specifications for maximum infusion pressure.
- Removal
  - Fix the guidewire using standard guidewire exchange techniques and carefully withdraw the Crossing Catheter

### Precautions

Store in a cool, dry place. Protect from direct sunlight and high temperature.

Use only appropriately sized ancillary device, as shown in the Specifications below.

Maximum Infusion Pressure: 360 psi (2482kpa)

Use the catheter prior to the "Use By" date specified on the package

The catheter should only be used by physicians qualified to perform percutaneous, vascular interventions.

Precautions to prevent or reduce clotting should be taken when any catheter is used in the vascular system. Use of systemic heparinization and heparinized saline solution should be considered.

Exercise care while handling the catheter during procedure to reduce the possibly of accidental damage, kinking or bending.

Manipulation of the catheter should only occur under fluoroscopy.

### Complications

Vascular catheterization and/or vascular intervention may result in complications including but not limited to:

- Catheter related: Adverse reaction to device materials, device malfunction/ fracture, infection, distal embolization, vessel thrombosis/ clot formation, vascular damage/ dissection requiring surgical repair, incompatibility with accessories

Complications usually associated with endovascular procedures and percutaneous access including:

- Unstable angina
- Hypo/hypertension
- Acute myocardial infarction
- Arrhythmia
- Death
- Hematoma
- Hemorrhage at puncture site
- False aneurysm formation
- Renal dysfunction

Warning: If the catheter is damaged, this product may cut into a blood vessel wall. Extreme caution needs to be taken when removing a damaged device. In the case of complications resulting from the removal of the entire system, stop immediately the procedure, and perform appropriate treatment at the discretion of the physician.

### Specifications

Model (Ref.)	Guide Compatibility	Guidewire Compatibility	Effective Length (cm)	Max Outer Diameter	Balloon Compatibility	Max Pressure psi(kpa)	Extendable Tip Diameter	Tip Extension
WGM35065CE	MIN 6F (2mm)	.035" (.89mm)	65 cm	.06" (1.5mm/4.6F)	MIN 5F (1.7mm)	360 (2482)	.050" (1.27mm)	0.20" (5mm)
WGM35090CE	MIN 6F (2mm)	.035" (.89mm)	90cm	.06" (1.5mm/4.6F)	MIN 5F (1.7mm)	360 (2482)	.050" (1.27mm)	0.20" (5mm)
WGM35135CE	MIN 6F (2mm)	.035" (.89mm)	135cm	.06" (1.5mm/4.6F)	MIN 5F (1.7mm)	360 (2482)	.050" (1.27mm)	0.20" (5mm)
WGM18090CE	MIN 5F (1.7mm)	.018" (.46mm)	90 cm	.05" (1.3mm/3.8F)	MIN 4F (1.3mm)	360 (2482)	.032" (.81mm)	0.20" (5mm)
WGM18135CE	MIN 5F (1.7mm)	.018" (.46mm)	135cm	.05" (1.3mm/3.8F)	MIN 4F (1.3mm)	360 (2482)	.032" (.81mm)	0.20" (5mm)
WGM18150CE	MIN 5F (1.7mm)	.018" (.46mm)	150cm	.05" (1.3mm/3.8F)	MIN 4F (1.3mm)	360 (2482)	.032" (.81mm)	0.20" (5mm)

### MANUFACTURER'S LIMITED WARRANTY

Manufacturer warrants that the Wingman Crossing Catheter is free from defects in material and workmanship when used by the stated Use By date and when package is unopened and undamaged immediately before use. Manufacturer's liability under this warranty is limited to replacement or refund of the purchase price of any defective Wingman Crossing Catheter. Damage to the Wingman Crossing Catheter caused by misuse, alteration, improper storage or handling, or any other failure to follow these Instructions for Use will void this limited warranty. THIS LIMITED WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. No person or entity, including any authorized representative or reseller of Manufacturer, has the authority to extend or expand this limited warranty and any purported attempt to do so will not be enforceable against Manufacturer.

Standard Symbol Legend										
LOT	Batch Code	Use By Date	Do not use if package is damaged	Do not reuse	Caution: See instructions for use	STERILIZED	Sterilized using irradiation	Consult instructions for use		
REF	Catalogue Number	Manufacturer	Sterilized using ethylene oxide	Do not resterilize	Rx ONLY Caution: Federal (USA) law restricts this device to sale by or on the order of a physician	EC REP	Authorized representative in the European Community	Non-pyrogenic		

Manufacturer:  
ReFlow Medical  
208 Avenida Fabricante #100  
San Clemente, CA 92672, USA

Authorized EU Rep  
MediMark Europe Sarl.  
11 rue Emile ZOLA B.P. 2332  
38033 Grenoble Cedex 2, France

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